Instruments for Endoscopic and Open Thoracic Spine Procedures







Thoracic spine surgery is often considered to be challenging. Teams specializing in this type of intervention treat discs, spinal trauma, deformations and even tumors.


There are, of course, hardly any similarities between the surgical treatment of a spinal trauma and the treatment of a spinal cord tumor. Irrespective of the different indications, some similarities nevertheless exist with regard to the surgical access and visualization for the thoracic region.

As part of a long and rich history of innovation in minimally invasive surgery, KARL STORZ has developed surgical access and visualization strategies exclusively for thoracic spine surgery. This product portfolio continues to grow thanks to co-operation with some of the world's most experienced surgeons and meanwhile covers the entire spectrum ranging from open surgery to keyhole interventions through to the percutaneous approach. Visualization strategies range from exoscopy to endoscopically assisted techniques through to interventions under full endoscopic visualization.

This brochure classifies all options according to product type (endoscopes, coagulation instruments, forceps, curettes, retractors etc.). This makes it easy and convenient for the surgeon to assemble a instrument set to suit individual needs or to upgrade existing equipment in order to extend the range of indications in the thoracic region and/or to make surgical access even less invasive.

Assembling an instrument set for endoscopic and open thoracic spine procedures

To assemble an instrument set, the user can follow the steps shown below and choose from the different instrument categories.



Knot Tier and Sponge Holder

Bone Punches


Forceps

Curettes, Dissectors, Hooks and Knives

Accessories

Holding Systems


Telescopes








26003 AE	ENDOCAMELEON® HOPKINS® Telescope, autoclavable, variable direction of view from 0° – 120°, adjustment knob for selecting the desired direction of view, fiber optic light transmission incorporated, color code: gold	10 mm	32 cm
----------	---	-------	-------





Note: Instruments with curved jaws or distally angled sheaths must be used with trocars with flexible cannulas!


For KARL STORZ trocars only the individual components are numbered.

The catalog number for the **entire assembled trocar,** as shown above in **bold** print, does not appear on the instrument.

40103 D

40103 X

Hooks, Specula, Dilation Rods






28090 SG Speculum	-	6 cm
--------------------------	---	------

26020 OM	Dilation Rod, for HEINKEL-SEMM Dilation Set 26020 HM and 26020 HP	-	-
----------	--	---	---

3 /

Retractors




28163 VR	ROSENTHAL Angled Retractor, with movable retractor blades, in three lengths, can be clamped to the holding arm with the KSLOCK		
2x 28163 VRM	Retractor Blade	12 cm	
2x 28163 VRK	Retractor Blade	14 cm	
2x 28163 VRG	Retractor Blade	18 cm	
2x 28272 HB	Articulated Stand, L-shaped	-	
2x 28172 HR	Rotation Socket, to clamp to the operating table	-	



SFENTACLE

28163 TRX	TENTACLE Spine Retractor, complete, with retractor frame, 2 retractor arms, X-ray transparent blades in 3 different lengths, light carrier with corresponding tools	
Complete set, cons	isting of:	
28163 TRF	Retractor Frame	-
28163 TRA	Retractor Arm	-
2x 28163 TRD	Drive Screw	-
2x 28163 TRB	Retractor Arm, extra	-
4x 28163 TBB	PEEK Retractor Blade	70 mm
4x 28163 TBS	PEEK Retractor Blade	50 mm
4x 28163 TBM	PEEK Retractor Blade	30 mm
6x 28163 TSB	Holder, for retractor blades	-
28163 TRL	Light Carrier	-
28163 TSL	Holder, for light carriers	-
28163 TSD	Screwdriver	-
28163 TSH	Positioning Instrument	-
39260 TRC	Sterilization Container	-





28163 BFP dist wid	ticulating Fan Retractor, tally bendable up to 30°, finger bracing dth 7.5 cm , autoclavable, for use with cars size 11	10 mm	34 cm
-----------------------	--	-------	-------

Unipolar Coagulation

				unipolar
	Art. No.	Description	Diameter	Length
				0
	33321 ML	CLICKLINE KELLY Dissecting and Grasping Forceps, rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock connector for cleaning, double action jaws, long		
	33321 UL	CLICKLINE REDDICK-OLSEN Dissecting and Grasping Forceps, rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock irrigation connector for cleaning, double action jaws, robust		
	33321 R	CLICKLINE Dissecting and Grasping Forceps, rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock irrigation connector for cleaning, double action jaws, right-angled	5 mm	36 cm
	34321 MS	CLICKLINE METZENBAUM Scissors , rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock connector for cleaning, double action jaws, curved		
<u></u>	34321 EK	CLICKLINE Hook Scissors, rotating, dismantling, with connector pin for unipolar coagulation, with irrigation connection for cleaning, single action jaws, tips of jaws not crossing		

10 11

unipolar

Garnery

KARL STORZ

Art. No.	Description	Diameter	Length
	Handle with Trumpet Valve, for suction or irrigation,		

Handle with Trumpet Valve, for suction or irrigation, 30804 autoclavable, coagulating suction tubes and 5 mm suction and irrigation tubes	-	-	
---	---	---	--

3077	UF Coagulation and Dissection Electrode, L-shaped, with connector pin for unipolar coagulation	5 mm	36 cm
------	--	------	-------


		•	Uternalny
26775 UE	Coagulation and Dissection Electrode, spatula-shaped, blunt, with connector pin for unipolar coagulation	5 mm	36 cm


Bipolar Coagulation

bipolar

Art. No.	Description	Diameter	Length
26276 HD	TAKE-APART® MANHES Bipolar Coagulation Forceps , width of jaws 3 mm		
26276 HK	TAKE-APART® SCHNEIDER Bipolar Grasping Forceps , with connector pin for bipolar coagulation, with movable inner sheath and non-retracting jaws, robust jaws	5 mm	33 cm
38361 ML	RoBi [®] KELLY Dissecting and Grasping Forceps, CLERMONT-FERRAND model, rotating, dismantling, with connector pin for bipolar coagulation, with LUER- Lock irrigation connector for cleaning, double action jaws, especially suitable for dissection	5 mm	36 cm

Suction and Irrigation Tubes







261	173 BN	Suction and Irrigation Tube, with lateral holes, with two-way stopcock for single hand control	5 mm	36 cm
-----	--------	--	------	-------


Clip Applicator and Needle Holder





30173 FAO	KOH Macro Needle Holder , dismantling, with LUER-Lock irrigation connector for cleaning, single action jaws, straight jaws, with tungsten carbide inserts, with ergonomic handle, axial, disengageable ratchet, ratchet position top, for use with suture material size 0/0 – 7/0	5 mm	33 cm
-----------	---	------	-------

Knot Tier and Sponge Holder




<u> <u>Burn</u></u>			.)	RL STOR2 Carrany
	32540 PT	Surgical Sponge Holder, self-retaining	10 mm	30 cm

Bone Punches


Art. No.	Description	Jaws	Round Sheath Diameter	Working Length
			So one	
	Demacur diamontling			


	28163 RKG	Rongeur, dismantling, graduated scale at the distal end	5 mm		
	28163 RGG	Rongeur, dismantling, graduated scale at the distal end	8 mm	10 mm	30 cm
L 7' 7' 7' 7'	28163 RKK	Rongeur, upturned 30°, dismantling, graduated scale at the distal end, for use with trocars size 11 mm	5 mm		




()	28163 PKS	KERRISON Bone Punch , dismantling, graduated scale at the distal end, 40° upbiting, not through- cutting	3 mm	10 mm	30 cm
P	28163 PKG	KERRISON Bone Punch , dismantling, graduated scale at the distal end, 90° upbiting, not through- cutting			

Forceps







Curettes

	Art. No.	Description	Diameter	Jaw Size	Working Length
			1010 85		KARI STORZ Cormany
2 2	28163 WC	Curette, sharp, spoon-shaped, oval		3.5 x 5 mm	
	28163 WD	Same		3.5 x 9 mm	
	28163 WE	Same		5 x 6.5 mm	
	28163 WF	Same		5.5 x 7.5 mm	
2 2 2	28163 WK	Same, round	10 mm	5 mm	30 cm
2 3	28163 WL	Same, oblong		2.5 x 6 mm	
21 31	28163 WM	Same, oblong		3.5 x 9 mm	
	28163 SBK	Curette, stirrup-shaped, for straightening end plates		-	



0	28163 GI	Spoon Curette, sharp		4 x 6 mm	
e	28163 GK	Same	5 mm	3 x 4 mm	36 cm
	28163 GL	Spoon Curette, sharp, spoon angled 45°		3 x 4 mm	

	Art. No.	Description	Diameter	Ring Size Diameter	Working Length
Q			.000.00		KARL STORZ Bermany

\bigcirc	28163 WN	Ring Curette, graduated, distal end slightly curved	10 mm	8 mm	30 cm
0	28163 WG	Ring Curette, graduated, distal end slightly curved	TOTIIII	6 mm	30 611

Dissectors, Hooks

Art. No.	Description	Diameter	Length
 		KARL STOP	22 Germany
28163 FDW	Dissector, distally curved 45°	2.6 mm	36 cm
 		XAR, STORE Growy 2008 10	0-5
26169 DO	Mikro-Messer n. BERCI , spitz, ausfahrbar,	5 mm	31 cm
			KARL STONZ On wany
28163 WH	Hook, 90°, blunt	10 mm	30 cm
			2
28163 GH	Palpation Hook, hook angled 90°, blunt	5 mm	36 cm
1 1 P		MALE STORY &	
28163 BT	Stumpfer Taststab, skaliert, zur Verwendung mit Trokaren der Größe 11 mm	5 mm	30 cm
			KARL STORZ Germany
28163 MS	Stößel, stumpf, gerieft, zum Positionieren eines Knochenspans	10 mm	30 cm

Raspatories, Chisels

Art. No. Description Diameter Length 28163 RR Raspatory 5 mm 36 cm 28163 KR Same 2 mm 2

	ALASTAZ Demo:
--	------------------

28163 MF	Raspatory, for dissection	7 mm	
28163 RP	Rasp, for coarsening the bone surface	10 mm	30 cm
28163 MS	Tappet, blunt, groved, for positioning the bone implant	10 mm	

1 2 3 4 5 6	-		MARK STOP2 Demon

	28163 MN	Chisel, graduated, double-edged grinding, straight	7 mm	30 cm
	28163 MM	Same	10 mm	30 CH


28163 BO	Stumpfes Osteotom, skaliert, mit Silikonhandgriff, zur Verwendung mit	10 mm	30 cm

Trokaren der Größe 11 cm

 28163 BE	Rippenköpfchen Elevatorium, skaliert, mit sichelförmiger Schneide, mit Silikonhandgriff, zur Verwendung mit Spekulum 28163 BSS oder Spekulum 28090 SG	10 mm	30 cm

Art. No.	Description	Length
	1 2 3 4 5 C KARL STORZ Germany	

Accessories





	Collar, for deflecting the mallet, for use with Chisel 28163 MM	-
--	--	---




Holding Systems



	riolang oystem, autolavable, with quor release coupling roleoor
	including:
28272 RGB	Rotation Socket, to clamp to the OR table
	Articulated Stand, L-shaped
	Clamping Jaw, clamping range 16.5 to 23 mm

Notes



It is recommended to check the suitability of the product for the intended procedure prior to use.

Consent to receive electronic information

Yes, I agree to receive future information by email at the following address:

Email	Name
Department / Practice	Street address
ZIP, Town	Signature
Lagree to my data being stored at KARL STORZ for this purpose. I can withdray	w my consent at any time and without giving reasons by emailing KARL STOR7

I agree to my data being stored at KARL STORZ for this purpose. I can withdraw my consent at any time and without giving reasons by emailing KARL STORZ at info@karlstorz.com. KARL STORZ will not make these data available to third parties.



WWW.KARLSTORZ.COM

KARL STORZ GmbH & Co. KG Mittelstraße 8, 78532 Tuttlingen, Germany Postbox 230, 78503 Tuttlingen, Germany Phone: +49 (0)7461 708-0 Fax: +49 (0)7461 708-105 E-Mail: info@karlstorz.com www.karlstorz.com

 KARL STORZ Endoscopy-America, Inc.

 2151 East Grand Avenue

 El Segundo, CA 90245-5017, USA

 Phone:
 +1 424 218-8100

 Phone toll free:
 800 421-0837 (US only)

 Fax:
 +1 424 218-8525

 Fax toll free:
 800 321-1304 (US only)

 E-Mail:
 info@ksea.com





