

Your Vision. Our Future

Rhino-Laryngo Videoscope

NEW OLYMPUS ENF TYPE VT2

Outstanding Observation Capability and Treatment Performance.

High quality imaging

Close focus for accuracy

2.0 mm channel for treatment

Multi-processor compatibility





High-resolution office therapy

High Image Quality, Close Focus Observation. from Suction to Biopsy, a Wide Range of Treatments are Possible.

NBI

imaging



High quality images

The ENF TYPE VT2 delivers some of the brightest, clearest images in a videoscope, making observation less strenuous. facilitating a detailed image for diagnosis.



NBI

Observation with Narrow Band Imaging™ (NBI) is available when used with the EVIS

close focus







2.0 mm viewing distance

The close-focus feature permits close-up observation of the tissue as near as 2 mm, enabling better observation of the mucosa

treatment



The suction mechanism lets you keep a clear view with easy push-button switching.

Accessories


20 different types of accessories are available, giving you a wide range of instruments to choose from.

Laser

HF

Compatible with both laser/electrosurgery treatments.

Accessory Lineup



Specifications



Optical system

Directions of view Depth of field

Insertion tube

Distal end outer diameter : Flexible tube outer diameter : Working length : 365 mm

Instrument channel

Minimum visible distance

Bending section

: Up 130°, Down 130°

Total length : 635 mm NBI : Available LD/YAG laser compatible / HF compatible

Compatible Video Systems







VISERA Video System



OTV-SI Compact OFFICE Video System

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP. OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
Posttach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany Postlach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 2009/ Hamburg, Germany OLYMPUS AMERICA INC.
3500 Corporate Parkway, P.O. Box 610, Center Valley, PA 18034-0610, U.S.A.
OLYMPUS SURGICAL & INDUSTRIAL AMERICA INC.
One Corporate Drive, Organopherus, New York 1968, U.S.A.

OLYMPUS LATIN AMERICA, INC. 5301 Blue Lagoon Drive, Suite 290 Miami, Florida 33126-2097, U.S.A.

OLYMPUS SINGAPORE PTE LTD 4918 River Valley Road #12-01/04, Valley Point Office Tower, Singapore 2465/3

OLYMPUS HONG KONG AND CHINA LIMITED

Room 1520-1527 Orean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS (BEJIING) SALES & SERVICE CO., LTD.
Room, 1202, NCI Tower, At2 Jianguomenwai Avenue Chaoyang District Beijing, 100022, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071 Moscow, Malaya Kajurkakan, 10 - Bill 1 1/2 Ducisi 1 1/2 Paris OLYMPUS AUSTRALIA PTY LTD
31 Gilby Boad Mount Waverley, VIC. 3149 Australia

www.olympus.com