

Monopolar
electrosurgical
pencils

Legend

Important notes


The authors and Erbe Elektromedizin GmbH have taken the greatest possible care while compiling the catalogue with electrosurgical accessories. Nevertheless, it is impossible to exclude the possibility that errors may be contained herein.

The information and recommendations given here may not be construed as constituting any basis for any claims against Erbe Elektromedizin GmbH. Should legal regulations stipulate liability, such liability shall be limited to intentional misconduct or gross negligence.

The data on the use of instruments is based on clinical experience. Medical science is subject to a process of permanent change based on research and clinical experience. Therefore the products are subject to alterations in design and technique. Because of production tolerances dimensions in the catalogue may marginally differ.

We confirm that none of our products contain latex, except the rubber straps, Art. No. 20592-009 and 20592-011.

Further product information



Please read carefully the corresponding notes on use before using any instrument.

Electrosurgery

Monopolar electrosurgical pencils

04	SINGLE USE electrosurgical pencils (ø 2.4 mm)	12	Electrosurgical pencils for T-Series
06	IES smoke evacuation	12	Electrosurgical pencils for epilation (ø 0.8 mm)
07	Premium electrosurgical pencils (ø 2.4 mm)	13	Special electrosurgical pencils for electrodes ø 4 mm
08	Premium electrosurgical pencils (ø 4 mm)	14	Electrosurgical pencils for arthroscopy (ø 2 mm)
10	ECO-LINE electrosurgical pencils (ø 2.4 mm)	15	Electrosurgical pencils for nerve test
11	Premium electrosurgical pencils with extension (ø 4 mm)		

Overview

Rate of applications

*PREMIUM
reusable (blue or black pencil)*



*ECO-LINE
reusable (grey pencil)*



*SINGLE USE
single use (white pencil)*



4 basic shapes

with 3 buttons



with 2 buttons



with rocker switch



without button



Connection

with integrated connecting cable

without connecting cable

Monopolar electro-surgical pencils

∅ 2,4

Single use electro-surgical pencils (∅ 2,4 mm)

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 3 m, with needle electrode



No. 20190-117

= 10

4.3 kVp

30190-259



STERILE R

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 3 m, with spatula electrode without protective holster



No. 20190-106

= 10

4.3 kVp

30190-212



STERILE R

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 3 m, with spatula electrode and protective holster



No. 20190-107

= 30

4.3 kVp

30190-212



STERILE R

Electrosurgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 3 m, with spatula electrode without protective holster



No. 20190-109

= 10

4.3 kVp

30190-212



STERILE R

Electrosurgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 3 m, with spatula electrode and protective holster



No. 20190-108

= 30

4.3 kVp

30190-212



STERILE R

Monopolar electro-surgical pencils

ø 2,4

Single use electro-surgical pencils (ø 2,4 mm)

Electrosurgical pencil without button, VIO, ICC, ACC, Standard

with connecting cable 3 m, with spatula electrode. The generator is activated by foot switch.



No. 20190-116

= 10

4.3 kVp

30190-259



STERILE R


Monopolar adapters for electro-surgical pencils with connecting plug International	Connection	Piece	No.
Monopolar adapter, VIO, ICC, ACC 	Standard (MO 9/5 mm)	1	20183-028
Monopolar adapter, T-Series 	Standard (MO 10/4 mm)	1	20183-029

Accessories for units and modules

IES smoke evacuation

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, \varnothing 10 / 22 mm, International

with connecting cable 3 m, with spatula electrode and integrated smoke evacuationshaft extension 100 mm.



No. 20321-028

= 25

4.5 kVp

30321-001



STERILE EO

Clip-on handle for Slim-Line electrosurgical pencils, tip short 12 mm

with evacuation tubing 3 m and connection \varnothing 22 mm, without electrosurgical pencil



No. 20321-007

= 25

30321-000



STERILE EO

Clip-on handle for Slim-Line electrosurgical pencils, tip long 100 mm

with evacuation tubing 3 m and connection \varnothing 22 mm, without electrosurgical pencil



No. 20321-020

= 25

30321-000



STERILE EO

Monopolar electrosurgical pencils

ø 2,4

Premium electrosurgical pencils (ø 2.4 mm)

Slim-Line electrosurgical pencil with 2 buttons, VIO, ICC, ACC, Standard

with connecting cable 4 m



No. 20190-105



= 1

4 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-067



= 1

6 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electrosurgical pencil with rocker switch, VIO, ICC, ACC, Standard

with connecting cable 4 m



No. 20190-104



= 1

4 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electrosurgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-076



= 1

6 kVp

max 95 °C



30190-179

max 138 °C

Monopolar electro surgical pencils

Ø 4

Premium electro surgical pencils (Ø 4 mm)

Slim-Line electro surgical pencil with 2 buttons, VIO, ICC, ACC, Standard

with connecting cable 4 m



No. 20190-065



= 1

4 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electro surgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-066



= 1

6 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electro surgical pencil with rocker switch, VIO, ICC, ACC, Standard

with connecting cable 4 m



No. 20190-074



= 1

4 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electro surgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-075



= 1

6 kVp

max 95 °C



30190-179

max 138 °C

Electro surgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-062



= 1

4.3 kVp

max 95 °C



30190-179

max 138 °C

VIO ReMode electro surgical pencil with 3 buttons, only for VIO 300 D

with connecting cable 4 m and MF plug



No. 20190-115



= 1

5 kVp

max 95 °C



30190-257

max 138 °C

Monopolar electrosurgical pencils

ø 4

Premium electrosurgical pencils (ø 4 mm)

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, Standard / International

without connecting cable



No. 20190-045



= 1

4.3 kVp

max 95 °C



30190-179

max 138 °C

Monopolar connecting cables	Connection	Length	Piece	No.
Monopolar connecting cable, VIO, ICC, ACC, Standard	Standard (MO 9/5 mm) Standard (MO 9/5 mm)	4 m 5 m	1 1	20192-127 20192-094
Monopolar connecting cable, Martin Standard	Martin Standard	4 m	1	20192-099
Monopolar connecting cable, VIO, ICC, ACC, International	International	4 m	1	20192-110

Electrosurgical pencil without button, VIO, ICC, ACC, Standard

with connecting cable 4 m. The generator is activated by foot switch



No. 20190-047



= 1

4 kVp

max 95 °C



30190-179

max 138 °C

Monopolar electrosurgical pencils

ø 2,4

ECO-LINE electrosurgical pencils (ø 2.4 mm)

ECO-Line electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-111



= 1

5 kVp

max 95 °C



30190-255

max 138 °C

ECO-Line electrosurgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-112



= 1



5 kVp

max 95 °C



30190-255

max 138 °C

Monopolar adapters for electrosurgical pencils with connecting plug International	Connection	Piece	No.
Monopolar adapter, VIO, ICC, ACC 	Standard (MO 9/5 mm)	1	20183-028
Monopolar adapter, T-Series 	Standard (MO 10/4 mm)	1	20183-029

Monopolar electro-surgical pencils

ø 4

Premium electro-surgical pencils (ø 4 mm)

Slim-Line electro-surgical pencil with rocker switch, VIO, ICC, ACC, Standard

with connecting cable 4 m, with extension 70 mm



No. 20190-090



= 1

2.5 kVp

max 95 °C



30190-179

max 138 °C

Slim-Line electro-surgical pencil with rocker switch, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m, with extension 70 mm



No. 20190-093



= 1

2.5 kVp

max 95 °C



30190-179

max 138 °C

Monopolar electrosurgical pencils

Electrosurgical pencils for T-Series

Electrosurgical pencil with 2 buttons, T-Series, Standard

for electrodes \varnothing 4 mm, without connecting cable



No. 20190-048



= 1

4.3 kVp

max 95 °C



30190-180

max 138 °C

Monopolar connecting cables

Monopolar connecting cable, T-Series



Connection

T-series

Length

4 m
5 m

Piece

1
1

No.

20192-095
20192-128

\varnothing 0,8

Electrosurgical pencils for epilation (\varnothing 0.8 mm)

Epilation handle without button, T-Series, Standard

for electrodes \varnothing 0.8 mm, with connecting cable 2 m. The generator is activated by foot switch.



No. 20190-008



= 1

1 kVp

max 95 °C



30190-197

max 138 °C

Monopolar electro-surgical pencils

ø 4

Special pencils for electrodes ø 4 mm

Suction and coagulation electrode, suction tube bendable, length 125 mm, with suction control

suction ø 3 mm (14 G, Gauge), connection for electro-surgical units ø 4 mm,
tube coupling from ø 3 mm to ø 5 mm



No. 20191-136



= 1




2.5 kVp

max 95 °C



30191-688

max 134 °C

Monopolar connecting cables	Connection	Length	Piece	No.
Monopolar connecting cable, VIO, ICC, ACC, Standard 	Standard (MO 9/5 mm)	4 m	1	20192-097
Monopolar connecting cable, T-Series 	T-series	4 m	1	20192-098
Monopolar connecting cable, VIO, ICC, ACC, International ø 8 mm (Bovie Jack) 	International ø 8 mm (Bovie Jack)	4 m	1	20192-111

Monopolar electrosurgical pencils

∅ 2

Electrosurgical pencils for arthroscopy (∅ 2 mm)

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, Standard / International

without connecting cable



No. 20190-082



= 1



4 kVp

max 95 °C



30190-196

max 138 °C

Monopolar connecting cables	Connection	Length	Piece	No.
Monopolar connecting cable, VIO, ICC, ACC, Standard 	Standard (MO 9/5 mm)	4 m	1	20192-127 20192-094
	Standard (MO 9/5 mm)	5 m	1	
Monopolar connecting cable, VIO, ICC, ACC, International 	International	4 m	1	20192-110

Electrosurgical pencil with 2 buttons, VIO, ICC, ACC, non-Erbe units, International

with connecting cable 4 m



No. 20190-098



= 1

4 kVp

max 95 °C



30190-196

max 138 °C

Monopolar electrosurgical pencils

Electrosurgical pencils for nerve test

Electrosurgical pencil for nerve test with 1 button, for silver ball electrode, NT 2, ICC 350 ZMK

with connecting cable 4 m



No. 20190-148



= 1

0.25 kVp

max 95 °C



30190-262

max 138 °C

Silver ball electrode for nerve test, ø 1 mm, length 40 mm

Shape / Forming

Piece

No.



straight

1

21191-052



Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Germany

Phone +49 7071 755-0
Fax +49 7071 755-179
info@erbe-med.com
erbe-med.com